

## Reinheldt R 7000 ICP-MS/MS

**Tandem Triple Quadrupole Mass Spectrometry**

*Uncompromising Precision. Triple Quadrupole Power.*



The **Reinheldt R 7000** represents the pinnacle of analytical performance, engineered for laboratories that demand absolute accuracy in the most complex matrices. By leveraging advanced triple quadrupole technology, the R 7000 eliminates spectral interferences, delivering detection limits and stability that conventional single quadrupole systems simply cannot reach.

## Technical Excellence: The Reinheldt Advantage

### 1. Advanced Multi-Mode Collision/Reaction Cell (CRC)

The R 7000 features a high-efficiency hexapole CRC, allowing for the flexible use of kinetic energy discrimination (KED) and specific reaction gases (such as H<sup>2</sup>, O<sup>2</sup>, or NH<sup>3</sup>). This ensures the complete removal of polyatomic interferences, providing a "clean" signal for even the most challenging analytes like Arsenic, Selenium, and Iron.

### 2. Superior Dual-Quadrupole Architecture

- **Q1 (First Quadrupole):** Operates as a high-precision mass filter, allowing only the target precursor mass to enter the reaction cell.
- **Q2 (Collision/Reaction Cell):** Facilitates controlled ion-molecule reactions.
- **Q3 (Analytical Quadrupole):** Filters the resulting product ions, ensuring a breakthrough in signal-to-noise ratio.

### 3. Ultra-Stable Ion Source & Interface

Our high-frequency, self-excited RF generator maintains an exceptionally stable plasma, even when directly aspirating high-matrix samples or organic solvents. The optimized interface cone design reduces deposition, extending maintenance intervals and ensuring long-term mass stability.

## Key Benefits for Your Laboratory

- **Exceptional Sensitivity:** Achieve parts-per-quadrillion (ppq) detection limits across the periodic table.
- **Matrix Versatility:** Seamlessly switch between environmental water, clinical samples, high-purity chemicals, and food safety testing.
- **High-Throughput Workflow:** Integrated automation compatibility and rapid stabilization times maximize your daily sample capacity.
- **Reinheldt Clarity Software:** An intuitive interface that simplifies method development—from automated tuning to complex data reporting.

## Applications

Industry	Capabilities
<b>Environmental</b>	Ultra-trace analysis of heavy metals in seawater and soil extracts.
<b>Semiconductor</b>	Analysis of ppt-level impurities in high-purity process chemicals.
<b>Food &amp; Pharma</b>	Compliance with USP <232>/<233> and stringent food safety regulations.
<b>Geology/Mining</b>	Precise isotope ratio measurements and rare earth element (REE) analysis.

## Strategic Solutions for Trace Analysis

### Eliminating Mass Overlaps

In traditional ICP-MS, unknown ions enter the cell and create by-product interferences. The R 7000's **MS/MS mode** adds a quadrupole filter before the cell. This ensures that only target ions enter, completely eliminating "ghost peaks" from matrix elements like Zirconium or Tungsten during Arsenic analysis.

### Ultra-Trace Sulfur & Phosphorus

Analyzing Sulfur is notoriously difficult due to  $O^2$  interference. The R 7000 uses mass-transfer technology to move  $^{32}S^+$  to  $^{32}S^{16}O^+$  at mass 48. Because Q1 has already filtered out background ions like  $^{48}Ti^+$  and  $^{48}Ca^+$ , the result is a pure, sub-ppt detection limit.

## Contact Information

### Reinheldt GmbH

Feldstrasse 14, D-86922 Eresing, Germany

**Email:** sales@reinheldtgmbh.de

**Phone:** +49 (0) 8193 9906 108

**Fax:** +49 (0) 8193 9961 490

**Hours:** 09:00 – 17:00 (CET)

*Precision Engineering for Analytical Leaders*